FRA — Dispositif Médical : VIS EAGLE 4.6 MM / VIS ARTHROSTAB 4.6 MM

REF : VIS EAGLE 38.98.XX

REF : VIS ARTHROSTAB 38.21.XX

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Matériau : titane ISO 5832-3

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ENG — Medical Device : EAGLE 4.6 MM SCREW / ARTHROSTAB 4.6 SCREW

REF: EAGLE 4.6 SCREW 38.98.XX

REF: ARTHROSTAB 4.6 SCREW 38.21.XX

(XX lenght from 20 to 60 mm)

Material: Titanium ISO 5832-3

Indications for use

The use of these screws is mainly in 3 indications:

Talo- navicular arthrodesis / Mediopod arthrodesis (M1 / C1 Lapidus) / MPGO arthrodesis.

Contra indications, Warning and possible adverse effects

Do not use in case of allergy to one of the material components. General surgical contra indication: poor quality bone due to osteoporosis. . .

The implants should never be reused. It's essential that the surgeon is fully conversant with the prescribed surgical

technique and that due consideration is given to the information provided below.

Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

The correct handing of implant is paramount. Contouring or altering the device should be avoided. Implants must not be modified or adapted.

Premature loading or muscular activity must be avoided as this may be problematic.

The surgeon should select the appropriate size and dimensions of the implant to suit the physiological functions of the patient and the type of osteosynthesis to be treated.

The product should NOT be implanted in the human body in contact with other metallic materials or different chemical composition as this may cause an electrolyte reaction and corrosion may occur.

The patient should be instructed to report any unusual changes around the operated site to his physician (e.g.: sign of infection) In case of non respect of those indications and contra indications, possible adverse effects included infection, subclincial nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, bending, or fracture of the device and/or the bone.

The patient must state that he/she is fitted with a metallic implant during MRI examinations or scans. The professional in charge of the examination shall ascertain the safety of the implant and patient, along with the risk of artefacts and data loss.

SURGICAL TECHNICAL

(The specific ancillary equipment is required for fitting ORTHO CAPE implants)

1. Talo-navicular arthrodesis

Patient in supine position + Pneumatic tourniquet. Medial incision centered on talo-navicular. Disassembly of the talo-navicular joint. Ablation of the cartilage manually using a distractor and chisel, Cauchoux and gouge pliers.

Verification of the dash of the two prepared surfaces. Installation of a stabilizing plate 4 self-supporting screws. Then, with the aid of the scapie, 1 pin provided with the EAGLE 4.6 mm / ARTHROSTAB 4.6 mm kit is inserted distally proximally from inside to out and from top to bottom under the plate.

The intraosseous measurement of this spindle is then carried out by subtraction using the scoop. The spindle is then bored about 2/3 of the length into the bone area with the wick, taking care to cross the joint to be treated under scopic control.

The stepped wick is then introduced to the stop to prepare the area of the screw head. Lastly, the EAGLE 4.6 mm / ARTHROSTAB 4.6 mm screw corresponding to the measurement minus 8 mm is introduced.

This EAGLE 4.6 mm / ARTHROSTAB 4.6 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well on both sides of the joint to be fused and ensures the compression of its surfaces.

The spindle can then be removed.

Immobilization by postoperative 90 ° removable boot.

2. (M1 / C1 Lapidus)

Patient in supine position + Pneumatic tourniquet.

Medial incision centered on the first joint of Lisfranc M1-C1 Disassembly of the M1-C1 joint is carried out Ablation of the cartilage.

Repositioning and correction of metatarsus varus Hold by a spit.

Installation of a plate 4 locking screw stabilization. Then, with the aid of the scapie, 1 pin provided with the EAGLE 4.6 mm / ARTHROSTAB 4.6 mm kit, is inserted distally proximally from inside to out and from top to bottom under the plate.

The intraosseous measurement of this spindle is then carried out by subtraction using the scoop.

The spindle is then bored about 2/3 of the length into the bone area with the wick, taking care to cross the joint to be treated under scopic control.

The stepped wick is then introduced to the stop to prepare the area of the screw head.

Lastly, the EAGLE 4.6 mm screw / ARTHROSTAB 4.6 mm screw corresponding to the measurement minus 8 mm is introduced. This EAGLE 4.6 mm screw / ARTHROSTAB 4.6 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well on both sides of the joint to be fused and ensures the compression of its surfaces.

The spindle can then be removed. Protection by postoperative rigid shoe type Podalux (DJO) or boot 90 °.

3. MPGO arthrodesis

Patient in supine position + Pneumatic tourniquet. We realize a medial approach centered on the MPGO of about 30mm

Opening of the capsule, we practice a complete and including external artholysis of MPGO

Then we perform a resection of the two articular surfaces to the Gouge forceps congruent dash.

The blocking position is then defined according to the usual criteria and temporarily stabilize by 1 pin supplied with the EAGLE 4.6 mm / ARTHROSTAB 4.6 mm set introduced from proximal to distal under scopic control.

The intraosseous measurement of this spindle is then carried out by subtraction using the scoop. The spindle is then bored about 2/3 of the length into the bone area with the wick, taking care to cross the joint to be treated under scopic control. The stepped wick is then introduced to the stop to prepare the area of the screw head.

Lastly, the EAGLE 4.6 mm screw / ARTHROSTAB 4.6 mm screw corresponding to the measurement minus 8 mm is introduced. This EAGLE 4.6 mm screw / ARTHROSTAB 4.6 mm screw is introduced by controlling its progression in scopie, so that the 2 threads of the screw is well on both sides of the joint to be fused and ensures the compression of its surfaces.

The proximal pin can then be removed. We can secure the assembly completed by 2 cross pins so as to eliminate any rotational movement satellite. Closing of the capsule after sesamoid refocusing with separate points.

Protection by postoperative rigid shoe type Podalux (DJO) or boot 90 °.

Do not use if the packaging is damaged

Device disposal: In the event of explantation, the implant should be given to a specialist department to ensure that it is disposed

according to strict environmental and health regulations. A product explanted for deficiency should be decontaminated and returned to the manufacturer.

Sterilization method: gamma irradiation

ORTHO CAPE 13 rue Mespoul,
31400 TOULOUSE - FRANCE.

Products should not be stored for a long time at temperatures below 10°C or above 33°C

The devices are not intended to be sterilized or reused. Resterilization or re-use can cause infections or product performance loss. Preserve from humidity to protect the readability of the label

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