FRA — Dispositif Médical : VIS M.I.S.S.® (Mini Invasive System Screw)

REF : 38.52.13 (13mm) – 38.52.15 (15mm) – 38.52.17 (17mm) – 38.52.19 (19mm) – 38.52.22 (22mm)

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ENG — Medical Device : M.I.S.S.® Screw (Mini Invasive System Screw)

REF : 38.52.13 (13mm) – 38.52.15 (15mm) – 38.52.17 (17mm) – 38.52.19 (19mm) – 38.52.22 (22mm)

Self-drilling screw, self-tapping, conical with compression with a new threading system allowing full use of the screw.

Material: titanium ISO 5832-3

Indications

M.I.S.S.® screws are designed to hold bone fragments in place during bone healing.

The M.I.S.S.® screw has been specially designed for minimally invasive surgery of the foot, especially for the fixation of osteotomies in chevron.

Contra indications, Warning and possible adverse effects

Do not use in case of allergy to one of the material components.

General surgical contra indication: poor quality bone due to osteoporosis...

The implants should never be reused. It's essential that the surgeon is fully conversant with the prescribed surgical technique and that due consideration is given to the information provides below. Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

The correct handing of implant is paramount. **Contouring or altering the device should be avoided.** Implants must not be modified or adapted. Premature loading or muscular activity must be avoided as this may be problematic.

The surgeon should select the appropriate size and dimensions of the implant to suit the physiological functions of the patient and the type of osteosynthesis to be treated.

The product should NOT be implanted in the human body in contact with other metallic materials or different chemical composition as this may cause an electrolyte reaction and corrosion may occur.

The patient should be instructed to report any unusual changes around the operated site to his physician (e.g.: sign of infection)

In case of non respect of those indications and contra indications, possible adverse effects included infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, bending, or fracture of the device and/or the bone.

The patient must state that he/she is fitted with a metallic implant during MRI examinations or scans. The professional in charge of the examination shall ascertain the safety of the implant and patient, along with the risk of artefacts and data loss.

SURGICAL TECHNICAL

(The specific ancillary equipment is required for fitting ORTHO CAPE implants)

- At the opening of the peelable blister preserve the adhesive labels of traceability (1 for renewal of order, 4 for patient's file).
- Set the rotation speed to about 500 rpm.
- Position the screw on the outer cortex pre-perforated by the square tip and start the engine while exerting pressure.
- Once the step engaged, the screw will flee forwards, and in compression, the screw will break.
- A safety of the embrittlement causes an early separation of the screw and its screw rod: Continue screwing with the manual screwdriver (14.30.12).
- The head of the inserted screw is flushing with the external cortex.

Do not use if the packaging is damaged.

Disposal: In the event of explantation, the implant should be given to a specialist department to ensure that it is disposed according to strict environmental and health regulations. A product explanted for deficiency should be decontaminated and returned to the manufacturer.

sterilization method: gamma irradiation

Ortho CAPE
172 impasse Lépine
82000 MONTAUBAN FRANCE

MD

Products should not be stored for a long time at temperatures below 10°C or above 33°C

The devices are not intended to be sterilized or reused. Restertization or re-use can cause infections or product performance loss. Preserve from humidity to protect the readability of the label

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