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ENG — MEDICAL DEVICE: VARISATION STAPLES

REF: 14.55.01 (90° - 10mm) –**14.55.02** (90° - 8mm)
14.55.03 (26° - 10mm) –**14.55.04** (26° - 8mm)

Indications for use

Tes staples are manufactured from inox 316L, ISO 5832-1. They are used for osteosynthesis and realignment of small osseus fragments, and phalangeal shortening osteotomy.

Contra indications, Warning and possible adverse effects

Do not use in case of allergy to one of the material components.

Not for use with nickel-allergic patients.

General surgical contra indication: poor quality bone due to osteoporosis...

The implants should never be reused. It's essential that the surgeon is fully conversant with the prescribed surgical technique and that due consideration is given to the information provides bellow.

Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

The correct handing of implant is paramount. **Contouring or altering the device should be avoided.** Implants must not be modified or adapted. Premature loading or muscular activity must be avoided as this may be problematic.

The surgeon should select the appropriate size and dimensions of the implant to suit the physiological functions of the patient and the type of osteosynthesis to be treated.

The product should NOT be implanted in the human body in contact with other metallic materials or different chemical composition as this may cause an electrolyte reaction and corrosion may occur.

The patient should be instructed to report any unusual changes around the operated site to his physician (e.g.: sign of infection)

In case of non respect of those indications and contra indications, possible adverse effects included infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, bending, or fracture of the device and/or the bone.

The patient must state that he/she is fitted with a metallic implant during MRI examinations or scans. The professional in charge of the examination shall ascertain the safety of the implant and patient, along with the risk of artefacts and data loss.

SURGICAL TECHNIQUE

(The specific ancillary equipment is required for fitting ORTHO CAPE implants)

- Make one distal hole in the distal fragment, directed forward and outward.

- Insert the staple with the staight staple driver(90°) or the oblique staple driver (26°)

- Finish with the impactor (90° or 26°)

Do not use if the packaging is damaged.

Device disposal: In the event of explantation, the implant should be given to a specialist department to ensure that it is disposed according to strict environmental and health regulations. A product explanted for deficiency should be decontaminated and returned to the manufacturer.

Sterilization method: gamma irradiation

ORTHO CAPE 13 rue Mespoul,
31400 TOULOUSE - FRANCE.

MD

Products should not be stored for a long time at temperatures below 10°C or above 33°C

The devices are not intended to be sterilized or reused. Resterilization or re-use can cause infections or product performance loss. Preserve from humidity to protect the readability of the label

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